


What Components Can I Produce With Each Bag System?

450 mL - CPDA-1 - Quad (with Diversion Arm)

Manufacturer ID:
BB*QCD456A2
ABRI ID: BG-WD 450 4 DNG


RBC expiration from date of
donation: 35 days



450 mL Optisol - Double (CPD) (with Sampling Arm)


Manufacturer ID:
BB*AGD456A3
ABRI ID: BG-WD 450 OPT
DN

RBC expiration from date of
donation: 42 days (with addition
of Optisol)



500 mL - Leukocyte Reduction Set (with CPD/Optisol)

Manufacturer ID:
BB*W6Q506A2
ABRI ID: BG-WD 500ml
RBC expiration from date of
donation:
35 days (without Optisol)
42 days (with addition of Optisol)



The terms double and quad refer to the number of bags with ports.