


# HemoTemp<sup>®</sup> Indicator

- An instant readout of the "core" temperature of the blood bag
- Helps you comply with A.A.B.B. temperature guidelines
- Makes it easy for staff to monitor blood bags prior to use
- Calibrated to indicate the current temperature of the blood therefore eliminating unnecessary waste caused by uncertainty about temperature



## CORE CORRELATED INDICATOR 4 REVERSIBLE INDICATORS



## COLOR-TEMPERATURE INTERPRETATION

- Temperature is in the lower number of the indicator range. (1°C, 4°C, 7°C, 10°C)
- Temperature is in the middle of the indicator range. (2°C, 5°C, 8°C, 11°C)
- Temperature is in the higher number of the indicator range. (3°C, 6°C, 9°C, 12°C)



**BIOSYNERGY, INC.** 1940 E. DEVON • Elk Grove Village, IL 60007  
TIME / TEMPERATURE / SOLUTIONS 800-255-5274 OR 847-956-0471 • www.biosynergyinc.com