



## Animal Blood Resources International

4989 Bird Drive  
Stockbridge, MI 49285  
Phone: 517-851-9855

### Laboratory Services Submission Form

Please use Adobe to **type** your responses. Do not use Google docs.

Please review submission requirements on next page before shipping samples.

Please **print** the completed form to include it with the sample submission. **Thank you!**

Date of Submission:

Account Number:

**Send results to:** Name:

Email Address:

Phone Number: (if we need to reach you):

**Bill to:** Company Name:

Address:

PO Number (if needed):

Invoicing Email Address:

### SAMPLE AND SHIPPING REQUIREMENTS:

In order to provide the best quality results, we request that you provide the sample amount requested. Samples should be shipped with blue ice packs. Wrap samples carefully for protection to prevent contact with ice packs. **SAMPLES SHOULD BE SENT BY OVERNIGHT COURIER.** Packages can only be received on Saturday by **SPECIAL** arrangement **ONLY**.

	Test Code	Collection Date	Client Sample Identification	Animal Name	Breed	D (Donor) or R (Recipient)
1						
2						
3						
4						
5						
6						
7						
8						
9						
10						

Additional comments:

☐ Check here if your company does not have an account with ABRI.

**Submitted by:**

## Laboratory Services Submission Requirements

To provide the best quality results, we request you provide the sample amount recommended.

### **Canine Tests** (Please call 517-851-9855 before sending crossmatch samples.)

#### Erythrocyte Testing

DEA	Minimum 1.5 mL EDTA whole blood <b>and</b> 2 mL separated serum
DEAS	Minimum 1.5 mL EDTA whole blood
DEA1S	Minimum 1.5 mL EDTA whole blood
DEA1	Minimum 1.5 mL EDTA whole blood <b>and</b> 2 mL separated serum
DEA4DEA7	Minimum 2 mL EDTA whole blood
DEA5	Minimum 2 mL EDTA whole blood
AbSCRN	Minimum 2 mL separated serum

#### Blood Donor Screening

K9PCRDONOR	Two tubes with minimum 1.5 mL EDTA whole blood <b>and</b> two tubes with 2 mL separated serum
K9PCRDONOR1	Two tubes with minimum 1.5 mL EDTA whole blood <b>and</b> two tubes with 2 mL separated serum
K9PCRSCRN	Minimum 1.5 mL EDTA whole blood <b>and</b> 2 mL separated serum
K9IFADONOR	Minimum 1.5 mL EDTA whole blood <b>and</b> 2 mL separated serum
K9IFADONOR1	Minimum 1.5 mL EDTA whole blood <b>and</b> 2 mL separated serum
K9DIFASCRN	Minimum 2 mL separated serum
K9CXM	Major - 1 mL EDTA DONOR whole blood <b>and</b> 2.5 mL RECIPIENT separated serum Minor - 1 mL EDTA RECIPIENT whole blood <b>and</b> 2.5 mL DONOR separated serum

### **Feline Tests** (Please call 517-851-9855 before sending crossmatch samples.)

#### Erythrocyte Testing

FEA	Minimum 1.5 mL EDTA whole blood <b>and</b> 1 mL separated serum
-----	---

#### Blood Donor Screening

FePCRDONOR	Two tubes with minimum 1.5 mL EDTA whole blood <b>and</b> two tubes with 1 mL separated serum
FePCR	Two tubes with minimum 1.5 mL EDTA whole blood <b>and</b> two tubes with 1 mL separated serum
FeCXM	Major - 1 mL EDTA DONOR whole blood and 2.5 mL RECIPIENT separated serum Minor - 1 mL EDTA RECIPIENT whole blood and 2.5 mL DONOR separated serum